



Passat

FRAMEWARMERS

Operating instructions

Passat [25] Passat [25E]

ORIGINAL-USER MANUAL

Contents:

1. General safety instructions
2. Startup and Operation
3. Intended Use
4. Restrictions
5. Foreseeable Misuse
6. Maintenance and Cleaning
7. Faults and Defects
8. Disposal
9. Technical Data
10. Imprint



Foreword and use of the device

Congratulations on your purchase of a Passat Framewarmer. You bought a device of the latest generation providing optimal heat distribution in the processing beaker as well as an elegant design. All devices of this series offer a very high level of operating safety and a long service life. For this unit to be of service to you for many years to come, please take the time to read this user manual carefully, especially before using the unit for the first time.

The Passat Framewarmer is specially designed for the requirements of opticians and hearing care professionals. It is used for heating of eyeglass frames or hearing aids for adjustments or the fitting of lenses in eyeglass frames.

1. General safety instructions

- Before using the device read the user manual thoroughly and always comply with the information provided
- Check if the rated voltage on the rating plate (base plate) is identical to your mains voltage
- Only connect the device to a properly installed shockproof socket
- Make sure to use the device only on a flat and dust-free surface
- Beaker and cover will heat during operation. Do not touch the hot surfaces. Also keep a safe distance to the air outlets at the side of the beaker and to the opening of the jet when using it
- There is a risk of body injury or property damage due to failure to comply with instructions given in this manual. Deyle GmbH will not assume liability for damages or defects due to failure to observe instructions contained in this user manual

2. Startup and operation

- Please observe the supplied installation instructions for the built-in type
Put the device into operation by using the main switch (1) – the switch will light during operation

Figure 1: main switch



- The device will start heating now. Temperature may vary depending on the period of use
For point heating please use the jet with opened cover. To remove the jet, please wait until it has completely cooled down. Otherwise there is the risk of burns
- This device generates heat. Device safety is only guaranteed if the device is used as intended (see Chapter 3)

3. Intended use

- The Passat Framewarmer must only be used in inner areas. The device must not be used in moist, flammable or explosive areas
- In case of overheating a built-in overheating-protection will cut the power supply to the heater. After cooling down the device will turn back on
- Do not drop screws or other smaller parts into the device while working. This could cause damage to the device
- Only use the jet with completely opened cover. If the cover touches the jet it could bend
- If you are not using the device for a longer period of time, turn it off and disconnect it from the mains
- Let the device completely cool down before you wrap it up or put it away
- During non-use store the device in a dry, clean place
- Deyle GmbH accepts no liability for any consequences from unintended use

4. Restrictions

- Only operate the device with rated voltage/frequency
- Mains fluctuations can cause temperature fluctuations of the device
- The device is intended for operation in upright position. Do not operate the device in horizontal or inclined positions
- This device can be operated by children from the age of 8 or older as well as by persons with reduced physical, sensory or mental abilities, lack of experience or knowledge, if they are under supervision or have been instructed about how to safely operate this device and have fully understood the resulting hazards.
- Children must not play with this device. Maintenance and Cleaning must not be carried out by children without supervision.

5. Foreseeable misuse

Please avoid the following points of misuse for a safe operation of your Passat Framewarmer. Observe the instructions given in this manual at all times!

Operation of a damaged device

- Do not operate the device if it is obviously damaged or has been dropped
- Do not operate the device with a damaged supply cable

Use of liquids on the device or in the surrounding areas

- Never use any liquids on the device, also with power switched off
- Never use any liquids in the surrounding area of the device
- Always avoid splash, hose or pressure water

Manipulation of the device

- Do not remove any parts
- Do not modify or convert the device
- Do not cover any ventilation openings
- Never completely cover the device!-→Fire Hazard!
- Never stick any objects through the protective grating!

6. Maintenance and cleaning

Always remove the plug before cleaning!!!

- The device requires no special maintenance
- In case of commercial use the device needs to be checked by an electrician periodically. This can also be carried out by your specialist dealer or Deyle GmbH. Please feel free to contact us (see Chapter 10/11) (Only applies to Germany – For other countries local statutory regulations and provisions must be observed)
- Always keep the device clean
- Only use a damp cloth to clean the device. The cloth must not drop!
- For stronger contamination you can use a cloth soaked with dishwashing liquid
- Do not use running water or pressure water for cleaning!
- Do not use solvent-containing or aggressive cleaners on the device. This could cause severe damage to the powder-coated surfaces as well as to the plastic parts
- Dust in the beaker may lead to odor development. This odor will stop after the dust was combusted. This is not a hazard or a defect

7. Faults and defects

- If the temperature of the device will strongly increase put the device out of operation right away
- In the event of faults or defects contact your local specialist dealer or our customer service
- Only qualified service personnel should remove the cover or attempt to service this device. We recommend a specialist dealer or our customer service for service and repair to achieve the best result. During warranty period the device must only be repaired by our customer service. Otherwise the warranty will become void

8. Disposal

- The device should be taken out of service and being properly disposed at the end of its service life
- We also take back the device and dispose it properly
- For the disposal the national statutory regulations and provisions in the country of installation must be complied with

9. Technical data

| | |
|---------------------------------------|---|
| Name: | Passat 25/25E |
| Dimensions incl. Beaker: | 190x190x265 mm (7,48x7,48x10,43 inch) |
| Built-in type: | 190x190x275 mm (7,48x7,48x10,83 inch) |
| Weight incl. Packaging: | 2,05/1,95 kg(4,52/4,3 lbs) |
| Rated Voltage/Frequency: | 230V/50Hz |
| Protection class: | I |
| Degree of protection: | IP20 |
| Heat output: | 220 W |
| Overheating protection: | Double Protection, in case of brief overheating a bimetal-fuse cuts out and turns the device back on after cooling down. In case of longer term overheating a thermal fuse will trip. It can only be used once and needs to be replaced |
| Emission sound pressure level: | LpA <70 dB(A) |
| Scope of delivery: | Passat Device, jet, manual |

11. Imprint

Deyle GmbH
Tulpenstraße 50
DE-71394 Kernen i. R.
Phone: +49 (0) 71 51 4 10 46
Fax: +49 (0) 71 51 4 68 06
E-Mail: mail@deyle-power.de
Internet: www.deyle-power.de
www.passat-heissluftgeraete.de/en/

Management: Dipl.-Ing. Dr. Florian Mehl, Dipl.-Ing. Dr. Horst Mehl
Eintragung beim Amtsgericht Stuttgart: HRB 262197
Ust.-Id.-Nr.: DE147319563

Passat Type 25/32 german patent pending

© 2018 Deyle GmbH All Rights reserved. This manual may not be duplicated without permission of Deyle GmbH

DEYLE

